



# A convenient and ideal monocular for use as a microscope or telescope

## NEITZ MONOCULAR 4×, 6×, 8×

These high quality monoculars fit easily in the palm of your hand and will fulfill your daily vision needs.

Their zoom focus range, from 20cm to infinity, is ideal for viewing scenery at a distance and for viewing small print or details.

By employing roof prisms, they are very compact providing you with an exceptionally clear view.

A close-up lens with a stand is also available for observing tiny objects or reading and writing for a person with low vision.

A kit consisting of 3 different essential monoculars and a close-up lens is recommended for professionals prescribing monoculars for people with special vision needs.

These monoculars are also ideal for use at school, work, travel, museum and hobbies.

### ● Specifications

Monocular type	4 × 12	6 × 16	8 × 20
Magnification			
at infinity	4×	6×	8×
at minimum focal distance	5.5×	7.6×	11×
with close-up lens	12×	20×	25×
Minimum focal distance (cm)	20	25	30
Exit pupil (mm)	3.0	2.7	2.5
Angular Field of view	12.5°	9.3°	7.0°
Field of view (m) at 1,000m	219	163	122
Weight (g)	55	65	80
Dimensions (mm)	31φ×58	31φ×70	31φ×98

Distributed by:

**NEITZ INSTRUMENTS CO., LTD.**

4F ICHIBANCHO COURT, 15-21 ICHIBANCHO, CHIYODA-KU, TOKYO 102-0082, JAPAN

PHONE: +81-3-3237-0552 FACSIMILE: +81-3-3237-0554

URL: <https://www.neitz.co.jp/en/>