

- ▶ "Fast Simple Accurate" Full Automatic Operation
- ▶ Reliable Assessment of CSF (Contrast Sensitivity Functionality) for Eye Professionals
- ▶ Reveal Functional Disability for e.g. Pre/Post-Cataract Surgery Diagnosis, LASIK, etc
- ▶ Brightness of Optotype :
Daylight Mode (200cd/sqm)
Twilight Mode (10cd/sqm)
plus Glare Testing Mode
(All LED Light Source)
- ▶ Three Contrast Levels : 100% 25% 10%
- ▶ Luminance Consistency and Aging Free Optotype
- ▶ Easy to Operate - Simple Operator Training
Anyone Can Learn
- ▶ Minimal Operator Fatigue, and Less Error
- ▶ Compact Design (Table Top Model as shown)
- ▶ Compatible with Laptop (via USB) : Easy to Save,
Downloadable and Printable of Job Summaries
and Normative Data
- ▶ Graphic Indication by Integrated Thermal Printer
- ▶ Reads More What VA Test Cannot



Contrast Sensitivity Acuity Tester



CAT-CP2