

Neitz Binocular Loupes BL series

Headband Type



NEITZ

Made by Neitz for Your Medical Field

NEITZ Headband Loupes

Excellent balance and comfort for all medical professionals.

- Hinge joints support fast & easy vertical adjustment.

- Multiple adjustment points for desirable fit.

- Flip-up Function: easy to see unmagnified view



- **Adjustable PD & working distance**
Revolvable barrels enable you to adjust the working distance for your perfect vision. User-oriented unique function, made by NEITZ.

- **LED Spot illuminators** with outstanding brightness and illumination range are available. (optional)



NSI-III



NSI-X



Yellow Filter
(Optional for NSI-III)

Specification

Model	Magnification	Working Distance	Field of View Diameter	Weight (incl. Headband)
Galilean Loupes				
BLS-1 HB	1.5x – 2.0x	550 mm – 350 mm (21.7 in – 13.8 in)	145 mm – 120 mm (5.7 in – 4.7 in)	251 g
BLS-2 HB	2.0x – 2.5x	550 mm – 350 mm (21.7 in – 13.8 in)	113 mm – 58 mm (4.5 in – 2.3 in)	246 g
BLS-3 HB	2.5x – 3.0x	550 mm – 350 mm (21.7 in – 13.8 in)	78 mm – 38 mm (3.1 in – 1.5 in)	240 g
BLD-3 HB	2.5x – 3.0x	350 mm – 250 mm (13.8 in – 9.8 in)	80 mm – 50 mm (3.2 in – 2.0 in)	235 g
Prismatic (Keplerian) Loupes				
BLP-4 HB	4.0x – 4.5x	550 mm – 350 mm (21.7 in – 13.8 in)	80 mm – 50 mm (3.2 in – 2.0 in)	258 g
BLP-6 HB	5.5x – 6.0x	550 mm – 350 mm (21.7 in – 13.8 in)	40 mm – 25 mm (1.6 in – 1.0 in)	261 g

Neitz Instruments Co., Ltd.

4F Ichibancho Court, 15-21, Ichibancho, Chiyoda-ku, Tokyo 102-0082, Japan
Tel : (+81) 3-3237-0552 Fax: (+81) 3-3237-0554 URL: <https://www.neitz.co.jp/en/>

